

DEKO 190 & DEKO 192

Bedpan Washer-Disinfectors

Protection for Patients, Safety for Staff



Pictured: DEKO 190 X2



Infection control you can trust

- Full removal of *C. difficile* (Independent study by HygCen Germany GmbH in 2025) and other high-risk pathogens
- A₀-value adjustable up to 13,000
- Flexible programs for different waste types, soiling levels, and outbreak scenarios

Hygienic by design

- Patented freshwater circulation
- Automatic disinfection of chamber, pipework, and spray jets in every cycle
- Optional "Non-Touch" automatic door for safer, hands-free use
- Minimizes manual handling and eliminates need for slop sinks

Smart and reliable with maximum ROI

- Built-in digital traceability, MDR + EN ISO 15883-1 & -3 certified
- Extensive life cycle; on average 16 years, up to 20+ years
- A low-waste, environmentally responsible alternative to macerators
- Reliable, energy-efficient, and user-friendly – designed for 24/7 use


Widest washable item range



Portable sanitary pans



Urine bottles and kidney dishes



Suction bottles, measuring jugs and various hollowware



Hospital bowls (patient wash bowls)



Raised toilet seats and toilet brushes / holders



Shoes



Similar products to the previous ones and used for similar purposes within patient care

DEKO 190 & DEKO 192

The Bedpan Washer-Disinfector That Sets the Standard

NEW!
DEKO 192
PASS-THROUGH



PREMIUM MDR- CERTIFIED PRODUCT

Trusted performance and full compliance
with EN ISO 15883-1 & -3

Unmatched washing power

- 600 L/min recirculation pump + 10 spray jets for exceptional cleaning
- High washing efficacy by patented freshwater recirculation
- Full removal of C.difficile

Large capacity/cycle:

- 2 bedpans + 2 covers and 4 urine bottles
- 2 large wash bowls and 4 urine bottles
- 8 urine bottles
- 4 pairs of shoes
- or any mixed load with ease



QUALITY FROM
FINLAND

Pictured: DEKO 190 X2

Loading examples

