

# DEKO Instrument Washer-Disinfectors

Engineered in Finland, trusted worldwide



DEKO 25



DEKO 260



DEKO 2000



DEKO D32

## Reliable, safe and compliant infection control


- Validated washing and thermal disinfection ensure maximum patient safety
- CE marked, ISO 15883 compliant
- Conforms to Medical Device Directive (MDR transition ongoing)
- Exceptional reliability and service life in critical environments, 24/7
- Proven performance: trusted in 40+ countries worldwide
- 24-month warranty for peace of mind

## Space- and cost-efficient performance

- Exceptional life-cycle value: low procurement and maintenance costs
- Space-efficient installation with no additional service access required
- Wide selection of models and sizes
- Modular racks enable fully customized solutions for a variety of instruments, including robotic surgery systems like da Vinci, CMR Versius, Medtronic Hugo
- Fast, validated cycles support efficient workflows

## Customizable with intelligent features

- Additional chemical pumps
- Pure water feed
- Sample drawing tap
- Conductivity measuring
- Network connectivity



Examples of special racks:  
MIS and Dental instruments

HOSPITALS • CSSDS • LABORATORIES • DENTAL & VETERINARY CLINICS • BEAUTY SALONS

**DEKO**  
MEDTECH

At your service:

dekomedtech.com  
sales@dekomedtech.com



QUALITY FROM  
FINLAND

Rev. 22/10/2025