

DEKO DC2200 GT Drying Cabinet

Efficient drying with increased capacity




QUALITY FROM
FINLAND


Why Choose DEKO DC2200 GT?

- **Fast and powerful drying**
Faster instrument turnaround
- **Energy-optimized process**
With heat exchanger
- **Fully adjustable programs**
1–180 min / 30–80°C
- **Up to 13 shelves / 26 DIN**
Maximize load capacity
- **Flexible access**
Load or unload at any stage;
program pauses/resumes automatically
- **Modular design**
Fits any application, environment, or workflow
- **Reliability you can trust**
Engineered for 24/7 use
- **Cost-efficient long lifecycle**
- **Simple plug-and-play installation**


Application

Time and temperature controls along with effective HEPA filtered airflow allow fast drying of reprocessable devices. Single door and pass-through options for optimizing workflows in CSSDs.


The DC2200GT Drying Cabinet is suitable for automatic drying of medical devices intended for re-use such as;



Surgical instruments and
MIS instruments



Instrument trays



Rigid endoscopes



Anaesthetic and
respiratory equipment



Transit containers



Hospital bowls and
other plasticware



Glassware



Hoses